

MX3 Ultrasound System

Transducer Family



C5-2m-e (Convex)

Applications: Abdominal, Vascular, Nerve, MSK, Small Parts, OB/GYN, Thoracic, Urology
Bandwidth: 1.2 MHz to 6.0 MHz
FOV: 68°
ExFOV: 80°
Radius: 50 mm
Needle guide: NGB-055, multi-angle, reusable



V10-4m (Endocavity)

Applications: OB/GYN, Urology
Bandwidth: 3.0 MHz to 12.0 MHz
FOV: 161°
Radius: 10 mm
Needle guide: NGB-004, NGB-045 single-angle, reusable



L12-3m-e (Linear)

Applications: Abdominal, Vascular, Nerve, MSK, Small Parts, Pediatric
Bandwidth: 3.0 MHz to 13.0 MHz
FOV (Max)
Non-extended: 3.79 cm
Extended: 10
Needle guide: NGB-007, multi-angle, reusable



L12-3Wm (Linear)

Applications: Abdominal, Vascular, Nerve, MSK, Small Parts, Pediatric
Bandwidth: 3.0 MHz to 13.0 MHz
FOV (Max):
Non-extended: 5.28 cm
Extended: 20
Needle guide: NGB-007, multi-angle, reusable



D6-2m (Volume Convex)

Applications: Abdominal, OB/GYN
Bandwidth: 1.8 MHz to 8.2 MHz
FOV (Max): 66 (B) x 86 (Sweep)
Radius: 40 mm
Needle guide: NA

Mindray North America

800 MacArthur Blvd., Mahwah, NJ 07430

Tel: 800.288.2121 Support: 877.913.9663 Fax: 800.926.4275 www.mindray.com

©2026 Mindray DS USA, Inc. Subject to change. Mindray® is a registered trademark of Shenzhen Mindray Bio-Medical Electronics Co. Ltd. 08/26 P/N: 0002-08-40908 Rev 2.0

mindray
healthcare within reach